



**Edit Deposition Time:**

- From main screen, click on "EDITORS" icon.
- **Focus on PM3 and select "Edit Step".** ←
- Select the appropriate "Category Filter", i.e. "Low-stress SiN Process(100% SiH4)". See list below for other options. ←
- Select the appropriate step from the "Directory of Process Steps". ←  
i.e. "SiN Dep 100C(100% SiH4)".
- Edit the "process Time" ← using the deposition rates on the paper posted on the right side of the system for your calculation and select "SAVE STEP". ← Select the "Return" icon. ←

**Category Filter**

**Step**


SiN process(2%SiH4-No-Ar)	{ SiN Deposit 50C-120W(No Ar) SiN Deposit 100C-120W(No Ar) SiN Deposit 250C-120W(NoAr)             }
"Low Stress No Ar"	
	{ SiN Deposition50C-50W(NoAr) SiN Deposition100C-50W(NoAr) SiN Deposition50C-50W(NoAr)             }
"Medium Stress No Ar"	
	{ SiN Deposition50C-5W(NoAr) SiN Deposition100C-5W(NoAr) SiN Deposition50C-5W(NoAr)             }
"High Stress No Ar"	

medium-stress SiN Process(100%SiH4)      SiN Dep 50C(100%SiH4)-50W  
SiN Dep 100C(100%SiH4)-50W  
SiN Dep 250C(100%SiH4)-50W

Low-stress SiN Process(100%SiH4)      SiN Dep 50C(100%SiH4)  
SiN Dep 100C(100%SiH4)  
SiN Dep 250C(100%SiH4)

SiO2 Process (100%SiH4)      "Low Dep Rate" {  
LDRSiO2 Dep 50C  
LDR SiO2 Dep 100C  
LDRSiO2 Dep 250C

"High Dep Rate" {  
HDRSiO2 Dep 50C  
HDR SiO2 Dep 100C  
HDRSiO2 Dep 250C



(Figure 2)

## SiNx and SiO2 with Argon Rate Table

### SiNx (Low Stress)

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	0	310	344
<u>Seconds/1000Å</u>	0	194	175

### SiNx (Medium Stress)

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	0	412	397
<u>Seconds/1000Å</u>	0	146	151

### SiO2 (LDR) - Low Dep Rate

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	375	382	383
<u>Seconds/1000Å</u>	160	157	157

### SiO2 (HDR) - High Dep Rate

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	1145	1137	1044
<u>Seconds/1000Å</u>	53	53	58

## SiNx and SiO2 with He Rate Table

### SiNx (Low Stress 2%SiH4 No Ar)

Flows: SiN-120W-100C(250C) (2%SiH4)

Category: SiN process (2%SiH4) Step: SiN deposition 100C( 50C, 250C)-120W

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	249	170	159
<u>Seconds/1000Å</u>	241	353	377

### SiNx (Medium Stress 2%SiH4 No Ar)

Flows: SiN-50W-100C(250C) (2%SiH4)

Category: SiN process (2%SiH4) Steps: SiN deposition 100C(50C,250C)-50W

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	289	219	212
<u>Seconds/1000Å</u>	207	274	283

### SiNx (High Stress 2%SiH4 No Ar)

Flows: SiN-5W-100C(250C) (2%SiH4)

Category: SiN process (2%SiH4) Steps: SiN deposition 100C(50C,250C)-5W

<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	300	237	209
<u>Seconds/1000Å</u>	200	253	287

### SiO2 (2%SiH4 with He)

Flows: SiO2-5W-100C(50C, 250C) (2%SiH4)








Category: SiO2 process (2% SiH4) Steps: SiO2 deposition 100C(50C, 250C)


<u>Substrate Temperature</u>	<u>50C</u>	<u>100C</u>	<u>250C</u>
<u>Dep Rate (Å/min)</u>	303	336	319
<u>Seconds/1000Å</u>	198	179	188

### Pre-Deposition Clean & Seasoning:

- From the "Lot Operations" screen select one of the clean & seasoning flows, and follow steps 6-7 of the "[Operation of Unaxis ICP-PECVD Deposition](#)" procedure using one of the following "Flows".
- To prepare for SiNx film at 250 °C, select "SF6 clean & SiN seasoning-250C".
- To prepare for SiNx film at 100 °C, select "SF6 clean & SiN seasoning-100C".
- To prepare for SiO2 film at 250 °C, select "SF6 clean & SiO2 seasoning-250C".
- To prepare for SiO2 film at 100 °C, select "SF6 clean SiNx & SiO2 seasoning-100C".
- Click on "Execute Sequence".

### PM3 Post-Dep Clean:


- Reference "Figure 4" and "Figure 5" on the next page for colored arrows within this procedure.
- From the "Edit Flow" screen  edit the "PM3 Post-Dep Clean" flow  and verify the Sequence "Post Dep Clean" is loaded into PM3 in the "Current Flow" list.  "Sequence" under "View Options" must be selected.
- From the "Edit Step" screen  select "Post Dep Clean" from the Category Filter list.  Select the step "SF6 Clean Post Dep"  and edit the "Process Time"  so the amount of material etched will equal your cumulative deposition thickness. SiNx etches at 20nm/min and SiO2 etches at 40nm/min.
- From the "Lot Operations" screen follow steps 6-7 of the "[Operation of Unaxis ICP-PECVD Deposition](#)" procedure using the flow "PM3 Post-Dep Clean".




(Figure 4)

### Editing a Flow:

- Navigate to the "Editors" screen.
- Select "Edit Flow" at the bottom of the screen.
- Select the "Flow" to be edited from the flow list.
- Click on "PM3" in the "Current Flow" field.
- Select the appropriate "Category Filter".
- Select the desired "Sequence" from the sequence list.
- If you are saving a custom flow make sure to rename the flow in the "flow name" field.
- "Save" the new flow if you are saving a standard flow like "PM3 Post-Dep Clean" then do not rename it.




(Figure 3)



(Figure 5)

## Aborting the run:

1. Wait until the wafer has loaded into the process chamber.
2. Select the "Stop Operations" icon. ←
3. Select the "STOP PM3 PROCESSING" button. ←
4. The "ALARMS" button will turn red, when it does select it. ←
5. Select "Next Step" then click on the "Recover" button.
6. Repeat steps 3-5 until the sequence reaches the "Final Pump Down" step.
7. Click on the "ALARMS" button to close the "ALARMS" screen. ←
8. Select the "MODULE STATUS" icon to return to the process status screen. ←



(Figure 6)