

# SURFSCAN 6200 for 8inch wafers

**IMPORTANT: Do not use it for wafers thicker than 1mm!**

1. **LOG IN** using the access code: **boss**  
(if the keyboard is not working use the virtual keyboard on the monitor)
2. Double click on the icon **SCAN**
3. Select the option **Recipe** (upper left corner) /**Load**
4. Load one of the Standard recipes for 8inch wafers  
**STD\_UCSB\_8INCH\_G2 -for (1.6 - 28) um**  
**STD\_UCSB\_8INCH\_G4 -for (0.16 - 1.6) um**  
**Do not make any changes in the recipes!**
5. Press OK, and recipe will be loaded  
(It will show the name of the recipe at the top of the monitor)
6. Select the option **Recipe**
7. Scroll down to option **Save as**
8. **Save as** - give a new name to the recipe (usually your name/company name etc.)
9. Select the option **Recipe/View/Modify**
10. **Modify** your recipe (wafer diameter, edge exclusion, gain, threshold, indexer and etc.)
11. On the same page select the option **Bin Split** and click on **Default Bin**  
This would update changes in your recipe.
12. **Save** the changes
13. **Loading wafers in the cassette:**  
(For 8" wafers use black cassette for loading the 8" wafers-left indexer)
14. Load the wafers into the cassette
15. Make sure wafer face is up, and major flat away from you
16. Place cassette in **the left indexer** (for 8")
17. Make sure cassette is sitting properly
18. Select on the monitor option **CASS** (Cassette)  
(Wait for the tool to read all wafers in the cassette.)
19. If you want to scan **ONLY ONE** wafer, select the wafer and click on **START** (upper bar)
20. Wait for the wafer to be picked up and scanned. Record results of the scan.
21. If you want to scan multiple wafers – select option **AUTO**, then **START**
22. When the scan is complete, click on the **HOME** option, select **RECIPE/exit**, and **LOG OFF**
23. Log in in the LOG BOOK and write details regarding your scan/or issues

**IMPORTANT: PLEASE LOG OFF after using the tool!**