

InP Etch using Unaxis PM1 at 200 C (2-17-2016)

Experimental:

- 1) PM1 Chamber O2 plasma clean for 30 minutes (Tony).
- 2) PM1 Chamber coating with a quarter of 2" InP wafer for 15 minutes.
- 3) Mounting a patterned InP sample onto a Si carrier with high-vacuum grease.
- 4) Sample Etching: 1.4mT, 125/800 W, Cl2/H2/Ar flow-rate=5/14/5 sccm, etch time=90 seconds.
- 5) Immediate after etching, soak the etched sample into DI water for a couple of minutes to remove chlorine residue.

Result:

Etch rate=1.11 $\mu\text{m}/\text{min}$; Etch selectivity (InP/SiO2) =13.4.

Figures 1 Etch profile.







